


OSNA[®] – fine-tune your radiation therapy with more differentiated information



Challenges in selecting the appropriate radiation therapy for your patient

Although guidelines recommend considering regional node irradiation (RNI) in early stage breast cancer patients with 1–3 positive sentinel lymph nodes (SLN), consensus is mixed among clinicians. Should one apply RNI in these patients? Which patients benefit? And which can be spared the treatment? There is also no clear recommendation as to the extent of the radiation required. Getting the dose right is currently a subjective decision.

More questions have been raised by clinical trials that investigated less invasive treatment approaches as a replacement for axillary lymph node dissection (ALND). Methodological flaws in all trials have led to remaining uncertainty about the appropriate radiation therapy (RT) for individual patients [1]. Moreover, in case of omission of ALND, decision-making is even more challenging as staging information regarding the nodal status is missing and has to be based solely on the result of the SLN. This becomes even more relevant as SLN analysis by histopathology is not standardised, examines only part of the lymph node and so provides only limited diagnostic information.




Abbreviations: BC: breast cancer; SLN: sentinel lymph node; ALND: axillary lymph node dissection; NACT: neoadjuvant chemotherapy; N: nodal; RNI: regional node irradiation.

Therefore, the question is: How do we address this information gap and is there a more standardised and comprehensive SLN assessment available that supports more reliable classification of patients?

OSNA[®] – more differentiated information to fine-tune treatment decisions

OSNA[®] – or One Step Nucleic Acid Amplification – is well established in routine and has proven its utility in many publications that show that the OSNA[®] result provides more differentiated diagnostic information than histopathology.



Overtreatment? No thanks!

The clinical trend is moving towards de-escalation of treatment in suitable patients to increase treatment efficacy and reduce unwanted side-effects. However, uncertainty remains about the appropriate loco-regional treatment for early breast cancer patients with no or low tumour burden in the SLNs. Should they be spared RNI and the related morbidity?

In this context, the importance of a standardised SLN assessment as well as the predictive [8] and prognostic [9] information as provided by OSNA® have become even more relevant.

The ongoing trials Optimal I and IIa aim to answer this question by exploring the correct axillary treatment for those patients in the adjuvant and neoadjuvant setting (Table 1, Figure 1). Both studies use OSNA® for lymph node analysis and were presented at the ESTRO 35 conference [10] and at the 'Assisi Think Tank Meeting' in which the importance of adherence and active participation to these trials was underlined during the meeting as well as in the resulting review [1].

Preliminary results of Optimal I presented at the AIS congress in Florence indicate similar outcome in both study arms.

Table 1

Trial	Type of study	Accrual time	Principal investigator, country	Cohort	Number of patients required	Main objective	Primary endpoint
Optimal I GIC-RAD-2014-0111	International, multicentre, prospective	April 2015 to December 2021	Manuel Ignacio Algara López, Spain	Early-stage BC, SLN assessed by OSNA® (250–15,000 copies/μL), no ALND	1,400	To demonstrate non-inferiority of incidental irradiation versus intentional irradiation to level I-III and supraclavicular nodes	5-year DFS
Optimal IIa GIC-RAD-2016-01	International, multicentre, prospective	January 2017 to January 2020	Manuel Ignacio Algara López, Spain	Early-stage BC, cN+ before NACT and ypNO after NACT (SLN assessed by OSNA®), no ALND	1,212	To demonstrate non-inferiority of irradiation to level I-II nodes versus irradiation to level I-III and supraclavicular nodes	5-year DFS

Abbreviations: ALND: Axillary lymph node dissection; DFS: Disease-free survival; NACT: neoadjuvant chemotherapy; SLN: Sentinel lymph node.




Figure 1 Schema of both Optimal I and II trials. Abbreviations mean BC: breast cancer; cNO: clinically negative lymph nodes; cN+: clinically positive lymph nodes; NACT: neoadjuvant chemotherapy; SLNB: sentinel lymph node biopsy; OSNA®: One Step Nucleic Acid Amplification; TTL: Total Tumour Load; ypNO: pathologically negative lymph nodes after NACT; RNI: regional node irradiation. * TTL range: 250–15,000 CK19 mRNA copies/μL.

References

- [1] **Aristei C, Kaidar-Person O, Arenas M, Coles C, Offersen BV et al.** (2016): *The 2016 Assisi Think Tank Meeting on breast cancer: white paper.* *Breast Cancer Res. Treat.* 160(2):211 – 221.
- [2] **Donker M, van Tienhoven G, Straver ME, Meijnen P, van de Velde CJH et al.** (2014): *Radiotherapy or surgery of the axilla after a positive sentinel node in breast cancer (EORTC 10981-22023 AMAROS): a randomised, multicentre, open-label, phase 3 non-inferiority trial.* *Lancet Oncol.* 15(12):1303 – 1310.
- [3] **Curigliano G, Burstein HJ, Winer EP, Gnant, Dubsy P et al.** (2017): *De-escalating and escalating treatments for early-stage breast cancer: the St. Gallen International Expert Consensus Conference on the Primary Therapy of Early Breast Cancer 2017.* *Ann. Oncol.* 28(8): 1700 – 1712.
- [4] **NCCN guidelines Breast Cancer 2017.** https://www.nccn.org/professionals/physician_gls/pdf/breast.pdf (Accessed 1 September, 2017).
- [5] **Giuliano AE, Ballman K, McCall L, Beitsch P, Whitworth PW et al.** (2016): *Locoregional Recurrence After Sentinel Lymph Node Dissection With or Without Axillary Dissection in Patients With Sentinel Lymph Node Metastases: Long-term Follow-up From the American College of Surgeons Oncology Group (Alliance) ACOSOG Z0011 Randomized Trial.* *Ann. Surg.* 264(3):413 – 420.
- [6] **Jagsi R, Chadha M, Moni J, Ballman K, Laurie F et al.** (2014): *Radiation field design in the ACOSOG Z0011 (Alliance) Trial.* *J Clin Oncol.* 2014, 32(32):3600 – 3606.
- [7] **Peg V, Espinosa-Bravo M, Vieites B, Vilardell F, Antúnez JR et al.** (2013): *Intraoperative molecular analysis of total tumor load in sentinel lymph node: a new predictor of axillary status in early breast cancer patients.* *Breast Cancer Res. Treat.* 139(1):87 – 93.
- [8] **Ohi Y, Umekita Y, Sagara Y, Rai Y, Yotsumoto D et al.** (2012): *Whole sentinel lymph node analysis by a molecular assay predicts axillary node status in breast cancer.* *Br. J. Cancer.* 107(8):1239 – 43.
- [9] **Peg V, Sansano I, Vieites B, Bernet L, Cano R et al.** (2017): *Role of total tumour load of sentinel lymph node on survival in early breast cancer patients.* *Breast.* 33:8 – 13.
- [10] **Algara López MI.** (2016): *The axilla – less surgery, more radiotherapy? In Kirby A (Chair), Regional nodal irradiation for breast cancer. Symposium conducted at the ESTRO 35 conference, Turin, Italy.*
- [11] **OPTimizing irradiation through molecular assessment of lymph node (OPTIMAL).** <https://clinicaltrials.gov/ct2/show/NCT02335957> (Accessed 1 September, 2017).
- [12] **Algara López MI, Sanz J et al.** (2017): *Irradiazione incidentale versus intenzionale negli stadi iniziali del cancro della mammella e con basso carico tumorale linfonodale. Risultati preliminari dall'OPTIMAL trial. Poster presented at the AIS conference, Florence, Italy.*

Authorised representative: **Sysmex Europe GmbH**

Bornbarch 1, 22848 Norderstedt, Germany · Phone +49 40 52726-0 · Fax +49 40 52726-100 · info@sysmex-europe.com · www.sysmex-europe.com

Manufacturer: **Sysmex Corporation**

1-5-1 Wakinohama-Kaigandori, Chuo-ku, Kobe 651-0073, Japan · Phone +81 78 265-0500 · Fax +81 78 265-0524 · www.sysmex.co.jp

You will find your local Sysmex representative's address under www.sysmex-europe.com/contacts